





Truewin Instruments for dental implants are made of high quality stainless steel.

Our line of instruments include, but is not limited to, the following products:

#### Bone Grafting Instruments

- Bone Curettes (SPRATT) – designed for harvesting bone
- Bone Grinder – used in conjunction with a mallet to grind bone into appropriate particle sizes
- Bone Basin – weighted bottom design prevents tipping or spilling while mixing bone or other bone grafting materials
- Bone Compactors – feature serrated tips for compacting bone

#### Scalpels & Biopsy Punches

- Scalpels – offers convenience of a disposable with razor sharp cutting edge for superior cutting
- Biopsy Punches – designed at a steeper angle to ensure a finely honed, sharper edge
- Bone Spreading Chisels
- Freer Chisels – double or single cut tips used for widening ridges during implant preparations
- Graduated Chisels – feature reference markings for measurement during a bone ridge augmentation and are available in straight and curved patterns

#### Osteotomes

- Concave Osteotomes – break the floor of the maxillary sinus through the implant cavity and are available in curved (off set) and straight (with or without stops)
- Convex Osteotomes – raise the floor of the maxillary sinus and are available in curved (off set) and straight patterns (with or without stops)
- Progressive Osteotomes – designed to prepare the implant site in the maxillary area where bone density is 3 or 4
- Bone Spreading Osteotomes – sharp or blunt configurations for progressive circumferential expansion of the implant site

#### Sinus Lift Instruments

- Sinus Lift Curettes and Elevators – lightweight, wide diameter stainless steel hand instruments lift the sinus membrane for bone augmentation procedures

#### Gauges & Calipers


- Gauges & Calipers – various measuring devices available to determine appropriate length, width, and distance for implant procedures

#### Collagen Products

- Collagen wound dressing product offering includes products.

#### Sterilization Trays

- Sterilization Trays – designed to protect and store instruments during sterilization, transportation and storage





- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

## IMPLANT GAGES



**TW6-01**  
Multifunctional Meter  
Length 90 mm,  
different sizes,  
8-10-11mm,  
corrosion-resistant

**TW6-02**  
Multifunctional Meter  
Length 90 mm,  
different sizes,  
5-13-15-17 mm,  
corrosion-resistant

**TW6-03**  
Measuring Probe  
Length 170 mm,  
(6-8-10-14-16 mm),  
corrosion-resistant

**TW6-04**  
Implant  
Depth Gauge

**TW6-05**  
Implant  
Depth Gauge  
& Key

**TW6-06**  
Implant  
Depth Gauge

**TW6-07**  
Implant  
Depth Gauge

## IMPLANT GAGES



**TW6-08**  
Double Ended Osteotomy Depth Gauge

**TW6-09**  
Double Ended Osteotomy Depth Gauge

**TW6-10**  
Double Ended Osteotomy Depth Gauge

**TW6-11**  
Depth & Gingival Gauges



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان



**TW6-12**  
Clamp, Osteometer  
Length about 155 mm,  
ridge width measurement,  
two-sided scale,  
measuring  
range 0 - 25 mm



**TW6-13**  
Measuring Loop  
straight length 90 mm,  
measuring range  
0 - 40 mm,



**TW6-14**  
Measuring Loop  
angled, length 90 mm,  
measuring range  
0 - 40 mm,



**TW6-15**  
Measuring loop, Castroviejo  
angled, length 175 mm,  
measuring range 0 - 40 mm,  
with window



**TW6-16**  
Boley Gauge  
Stainless Steel

TRUEWIN  
INDUSTRIES



**TW6-17**  
Measuring Loop  
angled, length  
175 mm,  
measuring range  
0 - 40 mm,



**TW6-18**  
Measuring Loop  
straight length  
175 mm,  
measuring range  
0 - 40 mm,



**TW6-19**  
Fig. 1  
Iwanson  
spring caliper  
For metal and  
porcelain.



**TW6-20**  
Fig. 2  
Iwanson  
spring caliper  
For wax.



**TW6-21**  
Bone Thickness Gauge  
40mm, 180mm




**TW6-22**  
Gauge  
for measuring  
distances  
(0 - 90 mm),



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

Hollow Handles 9.5mm Length 170mm  
BLACK TITANIUM COLOUR COATED INSTRUMENTS



TW6-23  
Fig. 1

TW6-24  
Fig. 2

TW6-25  
Fig. 3  
Extra-Flex

TW6-26  
Micro-Mini

TW6-27  
Mini 1

TW6-28  
Mini 4  
Extra-Flex

TW6-29  
Fig. 4  
Extra-Flex

TW6-30  
Fig. 5  
Flexi-Thin

TW6-31  
Fig. 6  
Flexi-Thin

TW6-32  
Mini 3  
Extra-Flex



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

Hollow Handles 9.5mm Length 170mm  
BLACK TITANIUM COLOUR COATED INSTRUMENTS



TW6-33  
AB1  
Boghosian

TW6-34  
AB2  
Boghosian

TW6-35  
W3

TW6-36  
8A

TW6-37  
4/5  
Gregg

TW6-38  
Interproximal  
Carver

TW6-39  
Interproximal  
Carver, Long

TW6-40  
A6

TW6-41  
Contouring  
Instrument



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

Hollow Handles 9.5mm Length 170mm  
**BLACK TITANIUM COLOUR COATED INSTRUMENTS**

**BURNISHERS**

**FREEDMAN BURNISHERS**



**TW6-42**  
21B

**TW6-43**  
27/29

**TW6-44**  
2  
Ladmore

**TW6-45**  
3  
Ladmore

**TW6-46**  
Freedman  
"Duckhead"

**TW6-47**  
Freedman  
Small Contact  
Forming


**TW6-48**  
Freedman  
Large Contact  
Forming

**TW6-49**  
Small/  
Medium Ball  
Burnisher



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

Hollow Handles 9.5mm Length 170mm  
**BLACK TITANIUM COLOUR COATED INSTRUMENTS**



**TW6-50**  
**A**  
 Cosmetic  
 Contouring

**TW6-51**  
**B**  
 Cosmetic  
 Contouring

**TW6-52**  
**C**  
 Cosmetic  
 Contouring

**TW6-53**  
**D**  
 Cosmetic  
 Contouring

**TW6-54**  
**E**  
 Cosmetic  
 Contouring

**TW6-55**  
**F**  
 Cosmetic  
 Contouring

**TW6-56**  
 Marginal Ridge  
 & Embrasure  
 Shaping Instrument


**TW6-57**  
 Occlusal  
 Anatomy  
 Instrument

**TW6-58**  
 Composite  
 Packing  
 Instrument



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان


Hollow Handles 9.5mm Length 170mm  
BLACK TITANIUM COLOUR COATED INSTRUMENTS





- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

**TITANIUM PLASMA COLOUR COATED  
SINUS LIFT INSTRUMENTS HOLLOW HANDLE 9.5MM LENGTH 170MM**



**TW6-69**  
Sinus Lift Elevator

**TW6-70**  
Sinus Lift Elevator

**TW6-71**  
Sinus Lift Elevator


**TW6-72**  
Sinus Lift Elevator

**TW6-73**  
Sinus Lift Elevator



- Инструменты для имплантологии
- Instrumentologie
- Implantologie
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

**TITANIUM PLASMA COLOUR COATED  
SINUS LIFT INSTRUMENTS HOLLOW HANDLE 9.5MM LENGTH 170MM**



**TW6-74**  
Sinus Lift Elevator

**TW6-75**  
Sinus Lift Elevator


**TW6-76**  
Sinus Lift Elevator

**TW6-77**  
Sinus Lift Elevator



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

**TITANIUM PLASMA COLOUR COATED  
IMPLANT SINUS LIFT INSTRUMENTS HOLLOW HANDLE 9.5MM LENGTH 170MM**



**TW6-78**  
Acute 6577

**TW6-79**  
Sinus1

**TW6-80**  
Sinus2

**TW6-81**  
Sinus3

**TW6-82**  
Sinus4

**TW6-83**  
Sinus5


**TW6-84**  
Sinus6

**TW6-85**  
Sinus7



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

**TITANIUM PLASMA COLOUR COATED  
SINUS LIFT INSTRUMENTS HOLLOW HANDLE 9.5MM LENGTH 170MM**



**TW6-86**  
Sinuslift 1

**TW6-87**  
Sinuslift 2

**TW6-88**  
Sinuslift 3

**TW6-89**  
Sinuslift 4

**TW6-90**  
Sinuslift 5

**TW6-91**  
Sinuslift 6

**TW6-92**  
Membrane  
3Positioning 1

**TW6-93**  
Condenser 1




- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

**TITANIUM PLASMA COLOUR COATED INSTRUMENTS**

**Sinus Lift Instrument,**  
Flexible Sinus Lift instruments  
reduce the risk of perforation of the  
"Schneider's membrane"

**Sinus Lift Elevator,**  
9.5mm Straight grip  
Length 175 mm, Sinus lift elevators, dull,  
for gentle peeling and lifting the  
Schneider membrane

**Sinus lift instrument**  
9.5mm Straight grip Length 180 mm,  
Buser (B1, B2, B3) raspatories,  
quality medical stainless steel,  
corrosion-resistant



TW6-94  
Fig. 1

TW6-95  
Fig. 2

TW6-96  
Fig. 3

TW6-97  
Fig. 1

TW6-98  
Fig. 2

TW6-99  
Fig. 3

TW6-100  
B01

TW6-101  
B02

TW6-102  
B03



- Инструменты для имплантологии
- Инструменти для имплантології
- Implantologie
- Implantologia
- Implantologie
- زراعة الأسنان

## TITANIUM PLASMA COLOUR COATED INSTRUMENTS

Sinus lift instrument 9.5mm Straight grip Sinus lift elevators, blunt or sharply curved (H01, H02, K03, K04, K78, K77)



**TW6-103**  
Sinus lift instrument  
H01

**TW6-104**  
Sinus lift instrument  
H02

**TW6-105**  
Sinus lift instrument  
K03

**TW6-106**  
Sinus lift instrument  
K04

**TW6-107**  
Sinus lift instrument  
K78

**TW6-108**  
Sinus lift instrument  
K77



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

**TITANIUM PLASMA COLOUR COATED INSTRUMENTS**

**SINUS LIFT INSTRUMENT**

Sinus Lift Elevators, 9.5mm Hollow Handle length 170mm Blunt Curved Set of 6 Pcs



TW6-109

TW6-110

TW6-111

TW6-112


TW6-113

TW6-114



- Инструменты для имплантологии
- Instrumentologie
- Implantologie
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

Hollow Handle 9.5mm Length 170mm  
PERIODONTAL SURGERY SURGICAL CHISELS  
KRAMER NEVINS SINUS LIFT





- Инструменты для имплантологии
- Инструменти для имплантології
- Implantologie
- Implantologia
- Implantologie
- زراعة الأسنان

**Hollow Handle 9.5mm Length 170mm  
IMPLANT SURGICAL CHISELS**



TW6-126 13K/TG



TW6-127 13K/TG Titanium Plasma Coated Point



TW6-128 13K/TG



TW6-129 13KL/R Titanium Plasma Coated Point




TW6-130 13KL 4/6mm



TW6-131 Pt2 Titanium Plasma Coated Point  
A/V Regular



TW6-132 Pt1 Titanium Plasma Coated Point  
A/V Regular



TW6-133 Pt3 Titanium Plasma Coated Point Angel S/E



TW6-134 Pt4 Titanium Plasma Coated Point STR S/E



TW6-135 Pt5 Titanium Plasma Coated Point TG CVD S/E



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

**TITANIUM PLASMA COLOUR COATED INSTRUMENTS**  
**Hollow Handle 10mm**  
**MIXING SPATULA WITH SHOVEL**



**TW6-136**  
 Bone Carrier  
 Titanium  
 185mm



**TW6-137**  
 Raspator, Prichard  
 Length 180 mm,  
 10 mm hollow handle,  
 double-ended  
 10.5 mm | 4.5 mm



**TW6-138**  
 Raspator, Prichard  
 Length 180 mm,  
 10 mm hollow handle,  
 double ended



**TW6-139**  
 Raspator, Prichard  
 Length 180 mm,  
 10 mm hollow handle,  
 double-ended  
 10.5 mm | 4.5 mm



**TW6-140**  
 Mixing Spatula  
 with shovel  
 170mm



**TW6-141**  
 Mixing Spatula  
 with shovel  
 185mm



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

**BONE GRAFTING INSTRUMENTS COMPACTOR /TAMPERS**

Product offering includes instruments for harvesting, mixing and compacting bone Provides clinicians with the proper instrumentation for bone grafting procedures Weighted bone basin Prevents tipping and spilling of bone and other materials during mixing

**Hollow Handle 9.5mm Length 170mm**



**TW6-142**  
Graft Packer,  
Concave Tips  
DE 4/6mm

**TW6-143**  
Bone Spoon/  
Packer DE 3mm

**TW6-144**  
Graft Packer,  
Serrated Tips  
DE 6/8mm

**TW6-145**  
Graft Packer,  
Serrated Tips  
DE 3/4mm

**TW6-146**  
Graft Packer,  
Smooth Tips  
DE 5/8mm

**TW6-147**  
Titanium graft packer/  
placement site depth  
gauge 3.3mm & 4.0mm  
Diameter tips Marked  
at 8, 11, 13, 15 & 18mm

**TW6-148**  
Graft Material Spoon  
Concave Tip For  
Larger Grafts



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

**BONE GRAFTING INSTRUMENTS**  
**TITANIUM PLASMA COLOUR COATED**  
**Hollow Handle 9.5mm**



**TW6-149**  
 Blade and pusher  
 Length 180 mm,  
 for straight or angled,  
 blade: 15 x 5 mm,  
 pusher: 4 mm diameter  
 Straight

**TW6-150**  
 Blade and pusher  
 Length 180 mm,  
 for straight or angled,  
 blade: 15 x 5 mm,  
 pusher: 4 mm diameter  
 angular

**TW6-151**  
 Spoon and pusher  
 Length 180 mm,  
 for angled spoon:  
 10 x 5 mm, pusher:  
 diameter 3 mm

**TW6-152**  
 Tunnelierungs  
 Instrument  
 length 175 mm,

**TW6-153**  
 Tunnelierungs  
 Instrument  
 length 175 mm,  
 different designs

**TW6-154**  
 Tamper  
 Length 170 mm,  
 for bone material,  
 double-ended,  
 diameter 3.2 mm  
 (2-4-6-8-10) and diameter  
 2.5mm (2-4-6-8-10)



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

Pluggers

Hollow Handle 9.5mm Length 160 mm, with titanium tips for mixing and pouring of autologous, heterologous bone substitute material, diameter (3/4 or 6/8 mm)

BONE GRAFTING INSTRUMENTS




TW6-155 Ø 4mm  
TW6-156 Ø 5mm  
Compactor  
150mm



TW6-157  
Sinus Window  
Bone Tamper  
4 mm



TW6-158  
Sinus Window  
Bone Tamper  
7 mm



TW6-159  
6/8mm



TW6-160  
3/4mm



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

## SPOON AND PUSHER (TRINOVO)

Hollow Handle 9.5mm Length 180mm



**TW6-161**  
10x5mm  
Spoon and  
Pusher  
with safety  
stop

**TW6-162**  
10x5mm  
Spoon and  
Pusher  
with safety  
stop

**TW6-163**  
10x5mm  
Spoon and  
Pusher  
with safety  
stop

**TW6-164**  
Bone Graft/Carrier  
Titanium  
2-4-6-8-10mm

**TW6-165**  
15x10mm  
Spoon and  
Pusher  
with safety  
stop

**TW6-166**  
10x5mm  
Spoon and  
Pusher  
with safety  
stop

**TW6-167**  
10x5mm  
Spoon and  
Pusher  
with safety  
stop



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

## PERIOSTEALS ELEVATOR 10mm Hollow Handle



**TW6-168**  
Raspator,  
Bone Scraper  
Length 165 mm,  
10 mm hollow  
handle,  
double-ended 6  
mm / 4 mm,

**TW6-169**  
Raspator,  
Raspatory  
Length 180 mm,  
10 mm hollow  
handle,  
double-ended  
3 mm | 3 mm,

**TW6-170**  
Raspator,  
Raspatory  
Length 180 mm,  
10 mm hollow  
handle,  
double-ended  
3 mm | 3 mm,

**TW6-171**  
Raspator,  
Papilla  
Length 185 mm,  
10 mm hollow  
handle,  
double-ended  
3 mm | 2 mm,

**TW6-172**  
Raspator,  
Bone Scraper  
Length 175 mm,  
10 mm hollow  
handle,  
double-ended,

**TW6-173**  
Raspator,  
Micro-Freer  
Length 180 mm,  
10 mm hollow  
double-ended  
3 mm | 4 mm

**TW6-174**  
Raspator,  
Freer  
Length 180 mm,  
10 mm hollow  
handle,  
double-ended  
5 mm | 5 mm,



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Інструменти для імплантології
  - Implantologia
  - زراعة الأسنان

## PERIOSTEALS ELEVATOR 10mm Hollow Handle



**TW6-175**  
Raspator,  
Freer  
10 mm hollow  
handle,  
high-quality  
medical  
double ended

**TW6-176**  
Raspator,  
Molt  
Length 180 mm,  
10 mm Hollow  
grip,  
double-ended  
6 mm

**TW6-177**  
Raspator,  
Micro Freer  
Length 185 mm,  
10 mm hollow  
handle,  
double-ended

**TW6-178**  
Raspator,  
Micro Freer  
Length 185 mm,  
10 mm hollow  
handle,  
double-ended

**TW6-179**  
Raspator, Micro  
Raspator /  
Chisel  
Length 170 mm,  
10 mm hollow  
handle,  
double-ended

**TW6-180**  
Raspator, Micro  
Raspator /  
Chisel  
Length 175 mm,  
10 mm hollow  
handle,  
double-ended

**TW6-181**  
Raspator,  
Micro Raspator /  
Chisel  
Length 175 mm,  
10 mm hollow handle,  
double-ended  
3.5mm | 3.5 mm



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

## PERIOSTEALS ELEVATOR 10mm Hollow Handle



**TW6-182**  
Raspatory Micro Molt  
Length 175 mm,  
corrosion-resistant

**TW6-183**  
Raspatory Chipping  
Workers:  
Chisel 24G and 7.5 mm,  
length about 175 mm,  
corrosion-resistant

**TW6-184**  
Raspatory Molt  
Work ends:  
5.0 mm and 4.0 mm,  
length about 175 mm,  
corrosion-resistant

**TW6-185**  
Rasp, Molt M9  
Length 175 mm,  
corrosion-resistant

**TW6-186**  
Raspator, Wautier  
Length 175 mm,  
corrosion-resistant


**TW6-187**  
Raspator, Benque  
Length 170 mm,  
workers  
5 mm | 7 mm,  
corrosion-resistant

**TW6-188**  
Raspator, Nordberg  
Length 180 mm,  
10 mm hollow handle,  
double-ended  
6 mm | 4 mm



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

## Curettes Hollow Handle 10mm Length 170mm





- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantoloia
  - زراعة الأسنان

## TISSUE MEMBRANE FORCEPS



**TW6-198**  
Tissue Graft Forcep  
Serrated Teeth  
150mm



**TW6-199**  
Self Locking  
Grooved  
150mm



**TW6-200**  
Tissue Graft Forcep  
Serrated Teeth  
150mm



**TW6-201**  
Titanium Metal  
Locking Plier



**TW6-202**  
TC Gripper Carbide  
Tipped Tissue Forceps  
Serrated Carbide Tips  
145mm



**TW6-203**  
Adson Carbide  
Tipped Tissue Forceps  
Serrated Carbide Tips  
120mm



**TW6-204**  
Micor Adson  
Tissue Forcep 120mm  
1x2 Teeth - 0.5mm Tips



**TW6-205**  
Adson 7x7  
Tissue Forceps  
7x7 Teeth 120mm



**TW6-206**  
Adson 7x7 Curved  
Tissue Forceps  
7x7 Teeth 120mm



**TW6-207**  
Suture Pliers  
Suture Forceps with  
1.6mm Ring Tip  
160mm



**TW6-208**  
Debakey Atraumatic  
Tissue Forceps Straight  
Serrated Teeth  
160mm



**TW6-209**  
Debakey Atraumatic  
Tissue Forceps Curved  
Serrated Teeth  
160mm



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

## TISSUE MEMBRANE FORCEPS



**TW6-210**  
Russian Long  
Tissue Forceps  
Serrated Tip  
200mm



**TW6-211**  
Russian  
Tissue Forceps  
Serrated Tip  
150mm



**TW6-212**  
Tissue Punch Forceps  
To Remove Tissue  
After Using Punch  
150mm



**TW6-213**  
narrow  
135mm



**TW6-214**  
105mm  
**Membrane tweezers KUEHNE**



**TW6-215**  
wide  
135mm



**TW6-216**  
Pin-Tweezers  
length 150mm  
with lock



**TW6-217**  
without lock  
150mm



**TW6-218**  
Haltepinzette with lock  
155mm, for pins, towards  
the tips, etc



**TW6-219**  
Carry Forceps  
150mm



**TW6-220**  
Carry Forceps  
150mm




**TW6-221**  
Carry Forceps  
150mm



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Інструменти для імплантології
  - Implantoloia
  - زراعة الأسنان


## OSTEOTOME CONVEXES / CONCAVE / CONIQUES

CONVEXES



**TW6-222**  
SET OF 5 PIECES  
Convexes

CONCAVES




**TW6-223**  
SET OF 5 PIECES  
Concaves

CONIQUES



**TW6-224**  
SET OF 5 PIECES  
Coniques



	CONVEXES STRAIGHT	<b>TW6-222</b>	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONCAVES STRAIGHT	<b>TW6-223</b>	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONIQUES STRAIGHT	<b>TW6-224</b>	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm



- Инструменты для имплантологии
- Implantologie
- Implantologie
- زراعة الأسنان
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان


## OSTEOTOME CONVEXES / CONCAVE / CONIQUES



**TW6-225**  
SET OF 5 PIECES  
Convexes Curved

**TW6-226**  
SET OF 5 PIECES  
Concaves Curved

**TW6-227**  
SET OF 5 PIECES  
Coniques Curved



	CONVEXES CURVED	<b>TW6-225</b>	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONCAVES CURVED	<b>TW6-226</b>	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONIQUES CURVED	<b>TW6-227</b>	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm



- Инструменты для имплантологии
- Implantologie
- Implantologie

- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

**TREPHINE**

Burs with depth markings are suited for precise and simple transplantation of cylindrical bone grafts as well as for all other bone-extraction tasks.



TRUWIN INDUSTRIAL



**TW6-241**  
Short Implant Drills  
Set of 3 Pcs

**TW6-242**  
Long Implant Drills  
Set of 4 Pcs





- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان

## BONE SCRAPER STRAIGHT & CURVED



**TW6-244**  
Bone Scraper  
Straight

**TW6-245**  
Bone Scraper  
Curved

## TISSUE PUNCH STRAIGHT & CURVED



Tissue Punch  
Straight

Tissue Punch  
Curved

**TW6-246**  
2.0mm

**TW6-247**  
4.0mm

**TW6-248**  
5.0mm

**TW6-249**  
6.0mm

**TW6-250**  
7.0mm

**TW6-251**  
8.0mm

**TW6-252**  
2.0mm

**TW6-253**  
4.0mm

**TW6-254**  
5.0mm

**TW6-255**  
6.0mm

**TW6-256**  
7.0mm

**TW6-257**  
8.0mm



**TW6-258**  
Bone Mill



**TW6-259**  
Basins For Bone



**TW6-260**  
Bone Crusher  
Metal With Titanium Crusher



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

## KEYS FOR IMPLANT SCREW



TW6-261



TW6-262



TW6-263



TW6-264

## TORQUE WRENCH & RATCHET



TW6-265 Star Torque Wrench 6.35mm Universal



TW6-266 Round Torque Wrench 4.0mm Universal



TW6-267 Hex Driver With 4.0mm Driver



TW6-268 Star Torque Wrench 6.35mm Universal



TW6-269 Hex Torque Wrench 4.00mm Universal



TW6-270 Star Torque Wrench 6.35mm Universal



TW6-271 Star Torque Wrench 6.35mm Universal



TW6-272 Star Torque Wrench 6.35mm Universal



TW6-273 Star Torque Wrench 6.35mm Universal



TW6-274 Star Torque Wrench 6.35mm Universal



TW6-275 Star Ratchet Wrench 6.35mm Universal



TW6-276 Star Ratchet Wrench 6.35mm Universal



TW6-277 Star Ratchet Wrench 6.35mm Universal



TW6-278 Star Ratchet Wrench 4.0mm Universal



TW6-279 Round Ratchet Wrench 4.0mm Universal



TW6-280 Round Ratchet Wrench 3.35mm Universal



TW6-281 Star Ratchet Wrench 4.0mm Universal



TW6-282 Star Ratchet Wrench 6.35mm Universal



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان



**TW6-285**  
Bone Crusher



**TW6-283**  
Bone Crusher  
Titanium Metal Cutter  
Set of 2 Pcs  
Bone Crusher + Mallet

mallet



**TW6-286**  
Bone Crusher  
Bone Grinding Grain  
1.5mm, 2mm, 2.5mm



**TW6-284**  
Bone Crusher  
Titanium Metal Cutter  
Set of 2 Pcs  
Bone Crusher + Mallet

mallet



**TW6-287**  
Bone Crusher  
Bone Grinding Grain  
1.5mm, 2mm, 2.5mm



- Инструменты для имплантологии
- Instrumenti для имплантології
- Implantologie
- Implantologia
- Implantologie
- زراعة الأسنان

## BONE EXPANDER SCREW & SURGICAL KITS



**TW6-288**  
 Bone Universal Kit  
 Drill: 2mm, 2.8mm, 3.2mm, 3.65mm,  
 4.2mm, 4.5mm, 4.8mm, 5.2mm  
 Trephines Bit:  
 3.0mm, 4.0mm, 5.0mm, 6.0mm,  
 7.0mm, 8.0mm, 9.0mm, 10mm.  
 Driver:  
 2.5mm, 2.5mm Hand Pcs  
 1.25mm, 1.25mm Hand Pcs  
 Short Driver 2.5mm  
 Long Driver 2.5mm  
 1.25mm, Short Hex Driver Hand Pcs  
 1.25mm, Long Hex Driver Hand Pcs  
 Torque Wrench 4.0mm



**TW6-289**  
 Bone Expander Kit  
 Expander: 2.6mm, 3.0mm,  
 3.4mm, 3.8mm, 4.3mm  
 Ratchet Driver  
 Parallel Pin  
 Pilot Drill  
 Drill 1.4mm  
 Drill Driver 4mm  
 Hand Driver 4mm  
 Ratchet 4.0mm



**TW6-290**  
 Bone Universal Kit  
 Drill: 2mm, 2.8mm, 3.2mm, 3.65mm,  
 4.2mm, 4.5mm, 4.8mm, 5.2mm  
 Trephines Bit:  
 3.0mm, 4.0mm, 5.0mm, 6.0mm,  
 7.0mm, 8.0mm, 9.0mm, 10mm.

Bone Compression:  
 2mm, 3.3mm, 2.8mm, 4.2mm, 3mm, 3.5mm  
 Driver:  
 2.5mm, 2.5mm Hand Pcs  
 1.25mm, 1.25mm Hand Pcs  
 Short Driver 2.5mm  
 Long Driver 2.5mm  
 1.25mm, Short Hex Driver Hand Pcs  
 1.25mm, Long Hex Driver Hand Pcs  
 Torque Wrench 6.35mm 1 Pcs



**TW6-291**  
 Conical Drills: 20 Pcs  
 Tissue Punches:  
 3.0mm, 3.50mm, 4.0mm, 4.5mm,  
 5.0mm, 6.0mm, 6.50mm  
 Bone Expander:  
 1.6/3.8mm, 1.9/3.0mm, 2.1/3.6mm,  
 2.6/3.8mm, 2.9/4.3mm  
 Drill Adapter: 1 Pcs  
 1.25mm, 2.50mm Hex Driver Hand Pcs

Pilot Drill: 1.6mm  
 Drill Bit: 1.4mm, 1.6mm, 1.8mm, 2.0mm, 2.2mm  
 Ratchet Driver: 4mm/6.35, 1.25/2.50  
 Rimer: 3.2/4.7mm, 3.7/4.7mm  
 Parallel Pin: 4.5mm, 2.5mm  
 Double step Drill: 2.4/2.8mm  
 Depth Gauge: 1 Pcs  
 Torque Wrench 6.35mm Universal: 1 Pcs



- Инструменты для имплантологии
  - Implantologie
  - Implantologie
- Инструменти для имплантології
  - Implantologia
  - زراعة الأسنان

## BONE INJECTORS, COLLECTORS & SUCTION TUBE



**TW6-292**  
 Ø 2,5mm  
 Straight / Bend  
 180mm  
 I.D 2.5mm  
 O.D 3.5

**TW6-293**  
 Ø 3,5mm  
 Straight / Bend  
 180mm  
 I.D 3.5mm  
 O.D 4.5

**TW6-294**  
 Ø 4,5mm  
 Straight / Bend  
 180mm  
 I.D 4.5mm  
 O.D 5.5

**TW6-295**  
 Ø 6,0mm  
 Straight / Bend  
 180mm  
 I.D 6.0mm  
 O.D 7.0

**TW6-296**  
 Ø 2,5mm  
 170mm  
 I.D 2.5mm  
 O.D 3.5

**TW6-297**  
 Ø 3,5mm  
 170mm  
 I.D 3.5mm  
 O.D 4.5

**TW6-298**  
 Ø 4,5mm  
 170mm  
 I.D 4.5mm  
 O.D 5.5

**TW6-299**  
 Ø 6,0mm  
 170mm  
 I.D 6.0mm  
 O.D 7.0



**TW6-300**  
 175mm  
 Ø 20mm

**TW6-301**  
 180mm  
 Ø 27mm

**TW6-302**  
 coarse  
 Ø 5mm

**TW6-303**  
 fine  
 Ø 9mm  
 Suction tubes  
 to collect  
 the extracted  
 bone material  
 180mm

**TW6-304**  
 Suction  
 Tubes  
 180mm

**TW6-305**  
 Suction  
 Tubes  
 180mm

**TW6-306**  
 Suction  
 Tubes  
 180mm

**TW6-307**  
 Fig. 6FR  
 Frazier  
 150mm - Ø 2 mm

**TW6-308**  
 Ø 3 mm  
 length 180mm  
**TW6-309**  
 Ø 4 mm  
 length 180mm  
**TW6-310**  
 Ø 5 mm  
 length 180mm



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Инструменти для имплантології
- Implantologia
- زراعة الأسنان



**TW6-311**  
Suction Cannula Kit  
with 3 Needles



Suction Pipe



Suction Needle  
1mm



Suction Needle  
1.5mm



Suction Needle  
2mm



**TW6-312**  
Height Gauge



**TW6-313**  
Brush Stainless Steel



**TW6-314**  
Amalgam Gun